

Audit of model-flagged anomalies

Each page shows one alert the model raised, in clinical context (6-second window, raw mV on standard ECG paper grid).

Total alerts in this PDF: 4

Records covered: 208

Centre beats hit by these alerts (cardiologist annotations):

V PVC	2
N Normal	1
F Fusion (V+N)	1


Legend on each page:

Black trace	the ECG (raw mV)
Pink shaded box	the 2-second window the model decided on
Green tick	annotated NORMAL beat (cardiologist label)
Magenta tick	annotated ABNORMAL beat (cardiologist label)
Letter under tick	MIT-BIH beat symbol

How to use this PDF:

1. Look at the pink window: that is what the model flagged.
2. Look at the magenta/green ticks INSIDE the pink window.
3. If there is a magenta tick inside or right at the edge, the model probably caught a real abnormality.
4. If everything inside the pink window is green, the alert is likely a false alarm.
5. Tick CORRECT FLAG / FALSE ALARM / NOT SURE on each page.


Alert #1/4 — record 208 at t = 00:32.49



Record 208 (MIT-BIH ID, NOT a patient identifier)
Time in record 32.49 s (sample 11697)
Model residual 0.0591
Alert threshold 0.0412
Mode fixed_online
Centre beat N - Normal
Wall-clock 2026-04-27T21:14:02.991265+00:00

Cardiologist audit: ☐ CORRECT FLAG ☐ FALSE ALARM ☐ NOT SURE
Notes: _____

Alert #2/4 — record 208 at t = 01:09.63




Record 208 (MIT-BIH ID, NOT a patient identifier)
Time in record 69.64 s (sample 25068)
Model residual 0.0681
Alert threshold 0.0412
Mode fixed_online
Centre beat F - Fusion (V+N)
Wall-clock 2026-04-27T21:14:40.810340+00:00

Cardiologist audit: ☐ CORRECT FLAG ☐ FALSE ALARM ☐ NOT SURE

Notes: _____

Alert #3/4 — record 208 at t = 01:17.86



Record 208 (MIT-BIH ID, NOT a patient identifier)
Time in record 77.87 s (sample 28033)
Model residual 0.0493
Alert threshold 0.0412
Mode fixed_online
Centre beat V - PVC
Wall-clock 2026-04-27T21:14:48.886299+00:00

Cardiologist audit: ☐ CORRECT FLAG ☐ FALSE ALARM ☐ NOT SURE
Notes: _____

Alert #4/4 — record 208 at t = 01:25.53

